



BLI-based Assessment of RBD-hACE2 Binding Interactions

Sample Preparation Data

Accession	Experimental ID	Construct ID
SARS_RBD	mut1	RBD-mut1-FC(IgG1)
KU182965	mut2	RBD-mut2-FC(IgG1)
EF065508	mut3	RBD-mut3-FC(IgG1)
EF065505	mut4	RBD-mut4-FC(IgG1)
EF065506	mut5	RBD-mut5-FC(IgG1)
EF065507	mut6	RBD-mut6-FC(IgG1)
MN611519	mut7	RBD-mut7-FC(IgG1)
MG916903	mut8	RBD-mut8-FC(IgG1)
JQ989266	mut9	RBD-mut9-FC(IgG1)
JQ989267	mut10	RBD-mut10-FC(IgG1)
JQ989268	mut11	RBD-mut11-FC(IgG1)
JQ989269	mut12	RBD-mut12-FC(IgG1)
KF294371	mut13	RBD-mut13-FC(IgG1)
AC_000192	mut14	RBD-mut14-FC(IgG1)
FJ647219	mut15	RBD-mut15-FC(IgG1)

Accession	Experimental ID	Construct ID
FJ647226	mut16	RBD-mut16-FC(IgG1)
FJ647220	mut17	RBD-mut17-FC(IgG1)
FJ647223	mut18	RBD-mut18-FC(IgG1)
MW532698	mut19	RBD-mut19-FC(IgG1)
MT121216	mut20	RBD-mut20-FC(IgG1)
OK017848	mut21	RBD-mut21-FC(IgG1)
OK017850	mut21	RBD-mut21-FC(IgG1)
OK017851	mut21	RBD-mut21-FC(IgG1)
OK017847	mut21	RBD-mut21-FC(IgG1)
OK017849	mut22	RBD-mut22-FC(IgG1)
OK017814	mut23	RBD-mut23-FC(IgG1)
OK017816	mut23	RBD-mut23-FC(IgG1)
OK017817	mut23	RBD-mut23-FC(IgG1)
OK017821	mut23	RBD-mut23-FC(IgG1)
OK017815	mut24	RBD-mut24-FC(IgG1)

Accession	Experimental ID	Construct ID
OK017820	mut25	RBD-mut25-FC(IgG1)
OK017825	mut26	RBD-mut26-FC(IgG1)
OK017826	mut27	RBD-mut27-FC(IgG1)
OK017844	mut28	RBD-mut28-FC(IgG1)
OK017852	mut29	RBD-mut29-FC(IgG1)
OK017856	mut30	RBD-mut30-FC(IgG1)
OK017853	mut30	RBD-mut30-FC(IgG1)
OK017855	mut30	RBD-mut30-FC(IgG1)
OK017854	mut31	RBD-mut31-FC(IgG1)
OK017857	mut32	RBD-mut32-FC(IgG1)
OK017797	mut33	RBD-mut33-FC(IgG1)
KT444582	mut34	RBD-mut34-FC(IgG1)
KY417150	mut35	RBD-mut35-FC(IgG1)
KF367457	mut36	RBD-mut36-FC(IgG1)
KY417151	mut37	RBD-mut37-FC(IgG1)

SARS_RBD and SARS2_RBD denote the reference RBD sequences of SARS-CoV and SARS-CoV-2 used for experimental comparison.

Sample Preparation Data

Accession	Experimental ID	Construct ID
MW681002	mut38	RBD-mut38-FC(IgG1)
MG772933	mut39	RBD-mut39-FC(IgG1)
AY278489	mut40	RBD-mut40-FC(IgG1)
AY572034	mut41	RBD-mut41-FC(IgG1)
AY572035	mut42	RBD-mut42-FC(IgG1)
AY686863	mut43	RBD-mut43-FC(IgG1)
AY686864	mut44	RBD-mut44-FC(IgG1)
AY572038	mut44	RBD-mut44-FC(IgG1)
DQ071615	mut45	RBD-mut45-FC(IgG1)
MG772934	mut46	RBD-mut46-FC(IgG1)
DQ022305	mut47	RBD-mut47-FC(IgG1)
DQ084200	mut47	RBD-mut47-FC(IgG1)
DQ084199	mut48	RBD-mut48-FC(IgG1)
GQ153539	mut49	RBD-mut49-FC(IgG1)
GQ153540	mut50	RBD-mut50-FC(IgG1)

Accession	Experimental ID	Construct ID
GQ153541	mut51	RBD-mut51-FC(IgG1)
GQ153544	mut52	RBD-mut52-FC(IgG1)
GQ153545	mut52	RBD-mut52-FC(IgG1)
GQ153546	mut52	RBD-mut52-FC(IgG1)
KP886808	mut53	RBD-mut53-FC(IgG1)
KP886809	mut54	RBD-mut54-FC(IgG1)
KF294457	mut55	RBD-mut55-FC(IgG1)
EU371561	mut56	RBD-mut56-FC(IgG1)
EU371562	mut57	RBD-mut57-FC(IgG1)
EU371563	mut58	RBD-mut58-FC(IgG1)
KJ473815	mut59	RBD-mut59-FC(IgG1)
MZ937001	mut60	RBD-mut60-FC(IgG1)
MZ937002	mut61	RBD-mut61-FC(IgG1)
MZ937003	mut62	RBD-mut62-FC(IgG1)
MZ937004	mut63	RBD-mut63-FC(IgG1)

Accession	Experimental ID	Construct ID
MZ937000	mut64	RBD-mut64-FC(IgG1)
KC881005	mut65	RBD-mut65-FC(IgG1)
EPI_ISL_410538	mut66	RBD-mut66-FC(IgG1)
EPI_ISL_410539	mut67	RBD-mut67-FC(IgG1)
EPI_ISL_410540	mut68	RBD-mut68-FC(IgG1)
EPI_ISL_410541	mut69	RBD-mut69-FC(IgG1)
EPI_ISL_410721	mut70	RBD-mut70-FC(IgG1)
EPI_ISL_471467	mut71	RBD-mut71-FC(IgG1)
EPI_ISL_471470	mut72	RBD-mut72-FC(IgG1)
C_AA085190	mut73	RBD-mut73-FC(IgG1)
C_AA085189	mut73	RBD-mut73-FC(IgG1)
C_AA085191	mut74	RBD-mut74-FC(IgG1)
C_AA085192	mut75	RBD-mut75-FC(IgG1)
C_AA085194	mut76	RBD-mut76-FC(IgG1)
SARS2_RBD	mut77	RBD-mut77-FC(IgG1)


SARS_RBD and SARS2_RBD denote the reference RBD sequences of SARS-CoV and SARS-CoV-2 used for experimental comparison.

SDS-PAGE

(Gel concentration:4~20%)

	Non-Reduced									Reduced								No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
Samples	IPI	1	2	3	4	5	6	M		1	2	3	4	5	6	IPI								
Load	1.0 μg	1.0 μg	1.0 μg	1.0 μg	1.0 μg	1.0 μg	1.0 μg	5 μL		1.0 μg	1.0 μg	1.0 μg	1.0 μg	1.0 μg	1.0 μg	1.0 μg								
<div> <div>kDa</div> <div>185</div> <div>115</div> <div>80</div> <div>65</div> <div>50</div> <div>30</div> <div>25</div> <div>15</div> <div>10</div> </div>																		4	JKY023	RBD-mut1-FC(IgG1)	P308333	Receptor binding domain	102.18	>85.0
																		5	JKY023	RBD-mut2-FC(IgG1)	P308334	Receptor binding domain	104.84	>85.0
																		6	JKY023	RBD-mut3-FC(IgG1)	P308335	Receptor binding domain	104.36	>85.0


SDS-PAGE (Gel concentration:4~20%)

	Non-Reduced									Reduced								No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
Samples	IPI	1	2	3	4	5	6	M		1	2	3	4	5	6	IPI								
Load	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	5 µL		1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg								
kDa 185 115 80 65 50 30 25 15 10																		1	JKY023	RBD-mut4-FC(IgG1)	P308336	Receptor binding domain	104.50	>85.0
																		2	JKY023	RBD-mut5-FC(IgG1)	P308337	Receptor binding domain	104.48	>85.0
																		3	JKY023	RBD-mut6-FC(IgG1)	P308338	Receptor binding domain	104.44	>85.0
																		4	JKY023	RBD-mut9-FC(IgG1)	P308341	Receptor binding domain	113.28	>85.0
																		5	JKY023	RBD-mut10-FC(IgG1)	P308342	Receptor binding domain	113.28	>85.0
																		6	JKY023	RBD-mut11-FC(IgG1)	P308343	Receptor binding domain	113.34	>85.0

SDS-PAGE (Gel concentration:4~20%)


	Non-Reduced								Reduced							No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
Samples	IPI	1	2	3	4	5	6	M	1	2	3	4	5	6	IPI							
Load	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	5 µL	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg							
<div><div>kDa</div><div>185</div><div>115</div><div>80</div><div>65</div><div>50</div><div>30</div><div>25</div><div>15</div><div>10</div></div>																5	JKY023	RBD-mut7-FC(IgG1)	P308339	Receptor binding domain	104.68	>85.0
																6	JKY023	RBD-mut8-FC(IgG1)	P308340	Receptor binding domain	113.46	>85.0

SDS-PAGE (Gel concentration:4~20%)


Non-Reduced									Reduced						
Samples	IPI	1	2	3	4	5	6	M	1	2	3	4	5	6	IPI
Load	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	5 µL	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg
kDa															
185															
115															
80															
65															
50															
30															
25															
15															
10															

No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
1	JKY023	RBD-mut12-FC(IgG1)	P308344	Receptor binding domain	113.28	>85.0
2	JKY023	RBD-mut13-FC(IgG1)	P308345	Receptor binding domain	120.62	>85.0
3	JKY023	RBD-mut14-FC(IgG1)	P308346	Receptor binding domain	123.68	>85.0
4	JKY023	RBD-mut15-FC(IgG1)	P308347	Receptor binding domain	123.68	>85.0
5	JKY023	RBD-mut17-FC(IgG1)	P308349	Receptor binding domain	123.68	>85.0

SDS-PAGE (Gel concentration:4~20%)




SDS-PAGE (Gel concentration:4~20%)

	Non-Reduced								Reduced							No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
Samples	IPI	1	2	3	4	5	6	M	1	2	3	4	5	6	IPI							
Load	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	5 µL	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg							
<div><div>kDa</div><div>185</div><div>115</div><div>80</div><div>65</div><div>50</div><div>30</div><div>25</div><div>15</div><div>10</div></div>																1	JKY023	RBD-mut18-FC(IgG1)	P308350	Receptor binding domain	120.50	>85.0

SDS-PAGE (Gel concentration:4~20%)

	Non-Reduced								Reduced							No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
Samples	IPI	1	2	3	4	5	6	M	1	2	3	4	5	6	IPI							
Load	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	5 µL	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg							
																4	JKY023	RBD-mut19-FC(IgG1)	P308351	Receptor binding domain	102.08	>95.0


SDS-PAGE (Gel concentration:4~20%)

	Non-Reduced								Reduced							No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
Samples	IPI	1	N/A	N/A	N/A	5	6	M	1	N/A	N/A	N/A	5	6	IPI							
Load	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	5 µL	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg							
<div><div>kDa</div><div><div>185</div><div>115</div><div>80</div><div>65</div><div>50</div><div>30</div><div>25</div><div>15</div><div>10</div></div></div> 	6	JKY023	RBD-mut20-FC(IgG1)	P308352	Receptor binding domain	102.32	>85.0															

SDS-PAGE (Gel concentration:4~20%)

	Non-Reduced									Reduced								No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
Samples	IPI	1	2	3	4	5	6	M		1	2	3	4	5	6	IPI								
Load	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	5 µL		1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg								
kDa 185 115 80 65 50 30 25 15 10																		1	JKY023	RBD-mut22-FC(IgG1)	P308354	Receptor binding domain	102.08	>85.0
																		2	JKY023	RBD-mut23-FC(IgG1)	P308355	Receptor binding domain	98.48	>85.0
																		3	JKY023	RBD-mut24-FC(IgG1)	P308356	Receptor binding domain	98.48	>85.0
																		4	JKY023	RBD-mut25-FC(IgG1)	P308357	Receptor binding domain	98.48	>85.0
																		5	JKY023	RBD-mut26-FC(IgG1)	P308358	Receptor binding domain	98.52	>85.0
																		6	JKY023	RBD-mut27-FC(IgG1)	P308359	Receptor binding domain	98.52	>85.0

SDS-PAGE (Gel concentration:4~20%)



No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
4	JKY023	RBD-mut28-FC(IgG1)	P308360	Receptor binding domain	98.40	>85.0
5	JKY023	RBD-mut31-FC(IgG1)	P308363	Receptor binding domain	101.58	>85.0
6	JKY023	RBD-mut32-FC(IgG1)	P308364	Receptor binding domain	101.58	>85.0

SDS-PAGE


(Gel concentration:4~20%)

	Non-Reduced									Reduced								No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
Samples	IPI	1	2	3	4	5	6	M		1	2	3	4	5	6	IPI								
Load	1.0 μg	1.0 μg	1.0 μg	1.0 μg	1.0 μg	1.0 μg	1.0 μg	5 μL		1.0 μg	1.0 μg	1.0 μg	1.0 μg	1.0 μg	1.0 μg	1.0 μg								
<div>kDa</div> <div>185</div> <div>115</div> <div>80</div> <div>65</div> <div>50</div> <div>30</div> <div>25</div> <div>15</div> <div>10</div>																		1	JKY023	RBD-mut36-FC(IgG1)	P308368	Receptor binding domain	102.16	>85.0
																		2	JKY023	RBD-mut35-FC(IgG1)	P308367	Receptor binding domain	102.16	>85.0
																		3	JKY023	RBD-mut34-FC(IgG1)	P308366	Receptor binding domain	102.16	>85.0
																		4	JKY023	RBD-mut33-FC(IgG1)	P308365	Receptor binding domain	98.22	>85.0
																		5	JKY023	RBD-mut37-FC(IgG1)	P308369	Receptor binding domain	102.08	>85.0
																		6	JKY023	RBD-mut38-FC(IgG1)	P308370	Receptor binding domain	98.58	>85.0

SDS-PAGE (Gel concentration:4~20%)

	Non-Reduced									Reduced								No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
Samples	IPI	1	2	3	4	5	6	M		1	2	3	4	5	6	IPI								
Load	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	5 µL		1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg								
kDa 185 115 80 65 50 30 25 15 10																		1	JKY023	RBD-mut39-FC(IgG1)	P308371	Receptor binding domain	98.50	>85.0
																		2	JKY023	RBD-mut40-FC(IgG1)	P308372	Receptor binding domain	102.46	>85.0
																		3	JKY023	RBD-mut41-FC(IgG1)	P308373	Receptor binding domain	102.00	>85.0
																		4	JKY023	RBD-mut42-FC(IgG1)	P308374	Receptor binding domain	102.00	>85.0
																		5	JKY023	ACE2-His	P308410	Receptor	69.91	>85.0

SDS-PAGE (Gel concentration:4~20%)

	Non-Reduced									Reduced								No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
Samples	IPI	1	2	3	4	5	6	M	1	2	3	4	5	6	IPI									
Load	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	5 µL	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg									
kDa 185 115 80 65 50 30 25 15 10																	3	JKY023	RBD-mut45-FC(IgG1)	P308377	Receptor binding domain	98.54	>85.0	
																	4	JKY023	RBD-mut46-FC(IgG1)	P308378	Receptor binding domain	98.58	>85.0	
																	5	JKY023	RBD-mut47-FC(IgG1)	P308379	Receptor binding domain	98.54	>85.0	
																	6	JKY023	RBD-mut48-FC(IgG1)	P308380	Receptor binding domain	98.54	>85.0	

SDS-PAGE


(Gel concentration:4~20%)

	Non-Reduced									Reduced								No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
Samples	IPI	1	2	3	N/A	5	6	M	1	2	3	N/A	5	6	IPI									
Load	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	5 µL	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg									
<div>kDa</div> <div>185</div> <div>115</div> <div>80</div> <div>65</div> <div>50</div> <div>30</div> <div>25</div> <div>15</div> <div>10</div>																1	JKY023	RBD-mut66-FC(IgG1)	P308398	Receptor binding domain	102.08	>85.0		
																2	JKY023	RBD-mut68-FC(IgG1)	P308400	Receptor binding domain	102.08	>85.0		
																3	JKY023	RBD-mut49-FC(IgG1)	P308381	Receptor binding domain	98.54	>85.0		
																4	JKY023	RBD-mut50-FC(IgG1)	P308382	Receptor binding domain	98.54	>85.0		
																5	JKY023	RBD-mut51-FC(IgG1)	P308383	Receptor binding domain	98.54	>85.0		
																6	JKY023	RBD-mut52-FC(IgG1)	P308384	Receptor binding domain	98.54	>85.0		

SDS-PAGE (Gel concentration:4~20%)

	Non-Reduced									Reduced								No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
Samples	IPI	1	2	3	4	5	N/A	M		1	2	3	4	5	N/A	IPI								
Load	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	5 µL		1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg								
kDa 185 115 80 65 50 30 25 15 10																		1	JKY023	RBD-mut53-FC(IgG1)	P308385	Receptor binding domain	98.28	>85.0
																		2	JKY023	RBD-mut54-FC(IgG1)	P308386	Receptor binding domain	98.28	>85.0
																		3	JKY023	RBD-mut55-FC(IgG1)	P308387	Receptor binding domain	98.60	>85.0
																		4	JKY023	RBD-mut56-FC(IgG1)	P308388	Receptor binding domain	102.18	>85.0
																		5	JKY023	RBD-mut57-FC(IgG1)	P308389	Receptor binding domain	102.18	>85.0
																		6	JKY023	RBD-mut58-FC(IgG1)	P308390	Receptor binding domain	102.18	>85.0

SDS-PAGE (Gel concentration:4~20%)

	Non-Reduced								Reduced							No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
Samples	IPI	1	N/A	N/A	N/A	5	6	M	1	N/A	N/A	N/A	5	6	IPI							
Load	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	5 µL	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg							
kDa																						
6	JKY023	RBD-mut59-FC(IgG1)	P308391	Receptor binding domain	98.48	>85.0																

SDS-PAGE


(Gel concentration:4~20%)

	Non-Reduced									Reduced								No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
Samples	IPI	1	2	3	4	5	6	M		1	2	3	4	5	6	IPI								
Load	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	5 µL		1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg								
<div>kDa</div> <div>185</div> <div>115</div> <div>80</div> <div>65</div> <div>50</div> <div>30</div> <div>25</div> <div>15</div> <div>10</div>																		1	JKY023	RBD-mut60-FC(IgG1)	P308392	Receptor binding domain	102.24	>85.0
																		2	JKY023	RBD-mut61-FC(IgG1)	P308393	Receptor binding domain	98.50	>85.0
																		3	JKY023	RBD-mut62-FC(IgG1)	P308394	Receptor binding domain	102.24	>85.0
																		4	JKY023	RBD-mut63-FC(IgG1)	P308395	Receptor binding domain	98.50	>85.0
																		5	JKY023	RBD-mut64-FC(IgG1)	P308396	Receptor binding domain	102.26	>85.0
																		6	JKY023	RBD-mut65-FC(IgG1)	P308397	Receptor binding domain	101.62	>85.0


SDS-PAGE (Gel concentration:4~20%)

	Non-Reduced								Reduced							No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
Samples	IPI	1	2	3	N/A	N/A	6	M	1	2	3	N/A	N/A	6	IPI							
Load	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	5 µL	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg							
kDa																1	JKY023	RBD-mut67-FC(IgG1)	P308399	Receptor binding domain	102.08	>85.0
185																						
115																						
80																						
65																						
50																						
30																						
25																						
15																						
10																						


SDS-PAGE (Gel concentration:4~20%)

	Non-Reduced									Reduced								No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
Samples	IPI	1	2	3	4	N/A	6	M	1	2	3	4	N/A	6	IPI									
Load	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	5 µL	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg									
<div>kDa</div> <div>185</div> <div>115</div> <div>80</div> <div>65</div> <div>50</div> <div>30</div> <div>25</div> <div>15</div> <div>10</div>																	1	JKY023	RBD-mut29-FC(IgG1)	P308361	Receptor binding domain	101.60	>85.0	
	2	JKY023	RBD-mut43-FC(IgG1)	P308375	Receptor binding domain	102.00	>85.0																	
	3	JKY023	RBD-mut44-FC(IgG1)	P308376	Receptor binding domain	102.02	>85.0																	
	4	JKY023	RBD-mut21-FC(IgG1)	P308353	Receptor binding domain	102.16	57.8																	
	6	JKY023	RBD-mut69-FC(IgG1)	P308401	Receptor binding domain	102.20	>85.0																	


SDS-PAGE (Gel concentration:4~20%)

	Non-Reduced									Reduced								No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
Samples	IPI	1	2	3	4	5	6	M		1	2	3	4	5	6	IPI								
Load	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	5 µL		1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg								
kDa 185 115 80 65 50 30 25 15 10																		1	JKY023	RBD-mut70-FC(IgG1)	P308402	Receptor binding domain	102.30	>85.0
																		2	JKY023	RBD-mut71-FC(IgG1)	P308403	Receptor binding domain	102.26	>85.0
																		3	JKY023	RBD-mut72-FC(IgG1)	P308404	Receptor binding domain	102.30	>85.0
																		4	JKY023	RBD-mut73-FC(IgG1)	P308405	Receptor binding domain	102.28	>85.0
																		5	JKY023	RBD-mut74-FC(IgG1)	P308406	Receptor binding domain	102.20	>85.0

SDS-PAGE (Gel concentration:4~20%)

	Non-Reduced									Reduced							No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
Samples	IPI	1	2	3	4	5	6	M	1	2	3	4	5	6	IPI								
Load	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	5 µL	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg	1.0 µg								
<div><div>kDa</div><div>185</div><div>115</div><div>80</div><div>65</div><div>50</div><div>30</div><div>25</div><div>15</div><div>10</div></div>																	1	JKY023	RBD-mut75-FC(IgG1)	P308407	Receptor binding domain	102.20	>85.0
	2	JKY023	RBD-mut76-FC(IgG1)	P308408	Receptor binding domain	102.24	>85.0																
	3	JKY023	RBD-mut77-FC(IgG1)	P308409	Receptor binding domain	102.34	>85.0																

SDS-PAGE (Gel concentration:4~20%)



No	Project No.	Sample	Protein No.	Protein Type	MW (kDa)	Purity (%)
5	JKY023	RBD- mut30-FC(IgG1)	P308362	Receptor binding domain	101.58	>85.0

Kinetics Result

Work flow :

1. Instrument: GATOR (ProbeLife)
2. Probe: ProA (LN 11062401)
3. Immobilization of ligands to probe: The antibody was diluted to 15 nM with Q buffer, The antibody coupling experiment is carried out .
4. Dynamic parameter setting: Adjust antigen to 400 nM with Q buffer, after diluting to 200, 100 nM, the combination time to 120 s, the dissociation time to 180 s, the default reaction temperature to 30°C, and then detection.
5. Result criteria: $R^2 > 0.95$, system adaptability passed, and data is authentic.

Number	Construct ID	Protein Number	ACE2-His (P308410)				
			KD (M)	ka(1/Ms)	kd(1/s)	R ²	R _{max} (nm)
1	RBD-mut1-FC(IgG1)	P308333	3.20E-07	1.80E+05	5.77E-02	0.9981	0.8364
2	RBD-mut2-FC(IgG1)	P308334	N/A	N/A	N/A	N/A	N/A
3	RBD-mut3-FC(IgG1)	P308335	N/A	N/A	N/A	N/A	N/A
4	RBD-mut4-FC(IgG1)	P308336	N/A	N/A	N/A	N/A	N/A
5	RBD-mut5-FC(IgG1)	P308337	N/A	N/A	N/A	N/A	N/A
6	RBD-mut6-FC(IgG1)	P308338	N/A	N/A	N/A	N/A	N/A
7	RBD-mut7-FC(IgG1)	P308339	N/A	N/A	N/A	N/A	N/A
8	RBD-mut8-FC(IgG1)	P308340	N/A	N/A	N/A	N/A	N/A
9	RBD-mut9-FC(IgG1)	P308341	N/A	N/A	N/A	N/A	N/A

KD: dissociation constant/association constant; **Ka:** association constant; **Kd:** dissociation constant; **R2:** Linear fitting constant; **Rmax:** The highest response value for curve fitting

Please note that the 'Number' in the table does not directly correspond to the 'mut' number in the Sample ID. For example, Number 30 corresponds to a Sample ID of RBD-mut31-FC(IgG1).

Kinetics Result

Number	Sample ID	Protein Number	ACE2-His (P308410)				
			KD (M)	ka(1/Ms)	kd(1/s)	R ²	R _{max} (nm)
10	RBD-mut10-FC(IgG1)	P308342	N/A	N/A	N/A	N/A	N/A
11	RBD-mut11-FC(IgG1)	P308343	N/A	N/A	N/A	N/A	N/A
12	RBD-mut12-FC(IgG1)	P308344	N/A	N/A	N/A	N/A	N/A
13	RBD-mut13-FC(IgG1)	P308345	N/A	N/A	N/A	N/A	N/A
14	RBD-mut14-FC(IgG1)	P308346	N/A	N/A	N/A	N/A	N/A
15	RBD-mut15-FC(IgG1)	P308347	N/A	N/A	N/A	N/A	N/A
16	RBD-mut16-FC(IgG1)	P308348	N/A	N/A	N/A	N/A	N/A
17	RBD-mut17-FC(IgG1)	P308349	N/A	N/A	N/A	N/A	N/A
18	RBD-mut18-FC(IgG1)	P308350	N/A	N/A	N/A	N/A	N/A
19	RBD-mut19-FC(IgG1)	P308351	6.54E-08	1.29E+05	8.43E-03	0.9979	1.0361
20	RBD-mut20-FC(IgG1)	P308352	6.31E-08	1.10E+05	6.93E-03	0.9959	1.4676
21	RBD-mut21-FC(IgG1)	P308353	N/A	N/A	N/A	N/A	N/A
22	RBD-mut22-FC(IgG1)	P308354	2.32E-07	2.64E+05	6.13E-02	0.9976	0.9049
23	RBD-mut23-FC(IgG1)	P308355	N/A	N/A	N/A	N/A	N/A

KD: dissociation constant/association constant; **Ka:** association constant; **Kd:** dissociation constant; **R2:** Linear fitting constant; **Rmax:** The highest response value for curve fitting

Please note that the 'Number' in the table does not directly correspond to the 'mut' number in the Sample ID. For example, Number 30 corresponds to a Sample ID of RBD-mut31-FC(IgG1).

Kinetics Result

Number	Sample ID	Protein number	ACE2-His (P308410)				
			KD (M)	ka(1/Ms)	kd(1/s)	R ²	R _{max} (nm)
24	RBD-mut24-FC(IgG1)	P308356	N/A	N/A	N/A	N/A	N/A
25	RBD-mut25-FC(IgG1)	P308357	N/A	N/A	N/A	N/A	N/A
26	RBD-mut26-FC(IgG1)	P308358	N/A	N/A	N/A	N/A	N/A
27	RBD-mut27-FC(IgG1)	P308359	N/A	N/A	N/A	N/A	N/A
28	RBD-mut28-FC(IgG1)	P308360	N/A	N/A	N/A	N/A	N/A
29	RBD-mut29-FC(IgG1)	P308361	2.82E-07	1.90E+05	5.37E-02	0.9942	1.0686
30	RBD-mut31-FC(IgG1)	P308363	1.28E-07	2.76E+05	3.53E-02	0.9942	0.867
31	RBD-mut32-FC(IgG1)	P308364	1.17E-07	2.58E+05	3.01E-02	0.9918	1.0486
32	RBD-mut33-FC(IgG1)	P308365	N/A	N/A	N/A	N/A	N/A
33	RBD-mut34-FC(IgG1)	P308366	1.43E-07	2.18E+05	3.12E-02	0.9946	1.2666
34	RBD-mut35-FC(IgG1)	P308367	1.32E-07	2.36E+05	3.10E-02	0.9925	1.1653
35	RBD-mut36-FC(IgG1)	P308368	1.17E-07	2.04E+05	2.39E-02	0.9881	1.0267
36	RBD-mut37-FC(IgG1)	P308369	2.42E-07	2.20E+05	5.32E-02	0.9944	0.926

KD: dissociation constant/association constant; **Ka:** association constant; **Kd:** dissociation constant; **R2:** Linear fitting constant; **Rmax:** The highest response value for curve fitting

Please note that the 'Number' in the table does not directly correspond to the 'mut' number in the Sample ID. For example, Number 30 corresponds to a Sample ID of RBD-mut31-FC(IgG1).

Kinetics Result

Number	Sample ID	Protein number	ACE2-His (P308410)				
			KD (M)	ka(1/Ms)	kd(1/s)	R ²	R _{max} (nm)
37	RBD-mut38-FC(IgG1)	P308370	N/A	N/A	N/A	N/A	N/A
38	RBD-mut39-FC(IgG1)	P308371	N/A	N/A	N/A	N/A	N/A
39	RBD-mut40-FC(IgG1)	P308372	2.71E-07	1.79E+05	4.86E-02	0.9976	1.0881
40	RBD-mut41-FC(IgG1)	P308373	N/A	N/A	N/A	N/A	N/A
41	RBD-mut42-FC(IgG1)	P308374	N/A	N/A	N/A	N/A	N/A
42	RBD-mut43-FC(IgG1)	P308375	N/A	N/A	N/A	N/A	N/A
43	RBD-mut44-FC(IgG1)	P308376	N/A	N/A	N/A	N/A	N/A
44	RBD-mut45-FC(IgG1)	P308377	N/A	N/A	N/A	N/A	N/A
45	RBD-mut46-FC(IgG1)	P308378	N/A	N/A	N/A	N/A	N/A
46	RBD-mut47-FC(IgG1)	P308379	N/A	N/A	N/A	N/A	N/A
47	RBD-mut48-FC(IgG1)	P308380	N/A	N/A	N/A	N/A	N/A
48	RBD-mut49-FC(IgG1)	P308381	N/A	N/A	N/A	N/A	N/A
49	RBD-mut50-FC(IgG1)	P308382	N/A	N/A	N/A	N/A	N/A
50	RBD-mut51-FC(IgG1)	P308383	N/A	N/A	N/A	N/A	N/A

KD: dissociation constant/association constant; Ka: association constant; Kd: dissociation constant; R2: Linear fitting constant; Rmax: The highest response value for curve fitting

Please note that the 'Number' in the table does not directly correspond to the 'mut' number in the Sample ID. For example, Number 30 corresponds to a Sample ID of RBD-mut31-FC(IgG1).

Kinetics Result

Number	Sample ID	Protein number	ACE2-His (P308410)				
			KD (M)	ka(1/Ms)	kd(1/s)	R ²	R _{max} (nm)
51	RBD-mut52-FC(IgG1)	P308384	N/A	N/A	N/A	N/A	N/A
52	RBD-mut53-FC(IgG1)	P308385	N/A	N/A	N/A	N/A	N/A
53	RBD-mut54-FC(IgG1)	P308386	N/A	N/A	N/A	N/A	N/A
54	RBD-mut55-FC(IgG1)	P308387	N/A	N/A	N/A	N/A	N/A
55	RBD-mut56-FC(IgG1)	P308388	2.43E-07	2.09E+05	5.08E-02	0.9951	1.0592
56	RBD-mut57-FC(IgG1)	P308389	2.37E-07	2.02E+05	4.79E-02	0.9955	1.1715
57	RBD-mut58-FC(IgG1)	P308390	2.16E-07	2.26E+05	4.90E-02	0.989	0.9693
58	RBD-mut59-FC(IgG1)	P308391	N/A	N/A	N/A	N/A	N/A
59	RBD-mut60-FC(IgG1)	P308392	2.43E-08	1.46E+05	3.54E-03	0.9862	1.6085
60	RBD-mut61-FC(IgG1)	P308393	N/A	N/A	N/A	N/A	N/A
61	RBD-mut62-FC(IgG1)	P308394	2.43E-08	2.01E+05	4.90E-03	0.9864	1.3462
62	RBD-mut63-FC(IgG1)	P308395	N/A	N/A	N/A	N/A	N/A
63	RBD-mut64-FC(IgG1)	P308396	1.24E-07	2.34E+05	2.89E-02	0.9897	1.3029
64	RBD-mut65-FC(IgG1)	P308397	5.46E-08	1.28E+05	7.01E-03	0.9921	1.4697
65	RBD-mut66-FC(IgG1)	P308398	5.39E-08	1.32E+05	7.11E-03	0.9929	1.2938

KD: dissociation constant/association constant; **Ka:** association constant; **Kd:** dissociation constant; **R2:** Linear fitting constant; **Rmax:** The highest response value for curve fitting

Please note that the 'Number' in the table does not directly correspond to the 'mut' number in the Sample ID. For example, Number 30 corresponds to a Sample ID of RBD-mut31-FC(IgG1).

Kinetics Result

Number	Sample ID	Protein number	ACE2-His (P308410)				
			KD (M)	ka(1/Ms)	kd(1/s)	R ²	R _{max} (nm)
66	RBD-mut67-FC(IgG1)	P308399	5.14E-08	1.24E+05	6.40E-03	0.9924	1.6548
67	RBD-mut68-FC(IgG1)	P308400	1.34E-07	5.74E+05	7.72E-02	0.9525	0.1651
68	RBD-mut69-FC(IgG1)	P308401	2.26E-08	1.76E+05	3.98E-03	0.9886	1.2357
69	RBD-mut70-FC(IgG1)	P308402	5.22E-08	1.29E+05	6.73E-03	0.9907	1.3488
70	RBD-mut71-FC(IgG1)	P308403	5.22E-08	1.49E+05	7.77E-03	0.9953	1.2501
71	RBD-mut72-FC(IgG1)	P308404	4.98E-08	1.53E+05	7.63E-03	0.994	1.2736
72	RBD-mut73-FC(IgG1)	P308405	1.22E-07	1.47E+05	1.79E-02	0.9691	0.2624
73	RBD-mut74-FC(IgG1)	P308406	N/A	N/A	N/A	N/A	N/A
74	RBD-mut75-FC(IgG1)	P308407	N/A	N/A	N/A	N/A	N/A
75	RBD-mut76-FC(IgG1)	P308408	9.29E-07	5.31E+04	4.93E-02	0.9518	0.3157
76	RBD-mut77-FC(IgG1)	P308409	1.11E-07	1.52E+05	1.69E-02	0.9989	1.3295
77	RBD-mut30-FC(IgG1)	P308362	1.52E-07	2.22E+05	3.37E-02	0.9966	0.7542

KD: dissociation constant/association constant; **Ka:** association constant; **Kd:** dissociation constant; **R2:** Linear fitting constant; **Rmax:** The highest response value for curve fitting

Please note that the 'Number' in the table does not directly correspond to the 'mut' number in the Sample ID. For example, Number 30 corresponds to a Sample ID of RBD-mut31-FC(IgG1).

N/A: very weak/no binding

Sensorgrams




Fig 1.RBD-mut1-FC(IgG1)




Fig 2. RBD-mut2-FC(IgG1)




Fig 3. RBD-mut3-FC(IgG1)




Fig 4. RBD-mut4-FC(IgG1)

Sensorgrams




Fig 5. RBD-mut5-FC(IgG1)




Fig 6. RBD-mut6-FC(IgG1)




Fig 7. RBD-mut7-FC(IgG1)




Fig 8. RBD-mut8-FC(IgG1)

Sensorgrams




Fig 9. RBD-mut9-FC(IgG1)




Fig 10. RBD-mut10-FC(IgG1)




Fig 11. RBD-mut11-FC(IgG1)




Fig 12. RBD-mut12-FC(IgG1)

Sensorgrams




Fig 13. RBD-mut13-FC(IgG1)




Fig 14. RBD-mut14-FC(IgG1)




Fig 15. RBD-mut15-FC(IgG1)




Fig 16. RBD-mut16-FC(IgG)

Sensorgrams




Fig 17. RBD-mut17-FC(IgG1)




Fig 18. RBD-mut18-FC(IgG1)




Fig 19. RBD-mut19-FC(IgG1)




Fig 20. RBD-mut20-FC(IgG1)

Sensorgrams




Fig 21. RBD-mut21-FC(IgG1)




Fig 22. RBD-mut22-FC(IgG1)




Fig 23. RBD-mut23-FC(IgG1)




Fig 24. RBD-mut24-FC(IgG1)

Sensorgrams




Fig 25. RBD-mut25-FC(IgG1)




Fig 26. RBD-mut26-FC(IgG1)




Fig 27. RBD-mut27-FC(IgG1)




Fig 28. RBD-mut28-FC(IgG1)

Sensorgrams




Fig 29. RBD-mut29-FC(IgG1)




Fig 30. RBD-mut31-FC(IgG1)




Fig 31. RBD-mut32-FC(IgG1)




Fig 32. RBD-mut33-FC(IgG1)

Sensorgrams




Fig 33. RBD-mut34-FC(IgG1)




Fig 34. RBD-mut35-FC(IgG1)




Fig 35. RBD-mut36-FC(IgG1)




Fig 36. RBD-mut37-FC(IgG1)

Sensorgrams




Fig 37. RBD-mut38-FC(IgG1)




Fig 38. RBD-mut39-FC(IgG1)




Fig 39. RBD-mut40-FC(IgG1)




Fig 40. RBD-mut41-FC(IgG1)

Sensorgrams




Fig 41. RBD-mut42-FC(IgG1)




Fig 42. RBD-mut43-FC(IgG1)




Fig 43. RBD-mut44-FC(IgG1)




Fig 44. RBD-mut45-FC(IgG1)

Sensorgrams




Fig 45. RBD-mut46-FC(IgG1)




Fig 47. RBD-mut48-FC(IgG1)




Fig 46. RBD-mut47-FC(IgG1)




Fig 48. RBD-mut49-FC(IgG1)

Sensorgrams




Fig 49. RBD-mut50-FC(IgG1)




Fig 50. RBD-mut51-FC(IgG1)




Fig 51. RBD-mut52-FC(IgG1)




Fig 52. RBD-mut53-FC(IgG1)

Sensorgrams




Fig 53. RBD-mut54-FC(IgG1)




Fig 54. RBD-mut55-FC(IgG1)




Fig 55. RBD-mut56-FC(IgG1)




Fig 56. RBD-mut57-FC(IgG1)

Sensorgrams




Fig 57. RBD-mut58-FC(IgG1)




Fig 58. RBD-mut59-FC(IgG1)




Fig 59. RBD-mut60-FC(IgG1)




Fig 60. RBD-mut61-FC(IgG1)

Sensorgrams




Fig 61. RBD-mut62-FC(IgG1)




Fig 62. RBD-mut63-FC(IgG1)




Fig 63. RBD-mut64-FC(IgG1)




Fig 64. RBD-mut65-FC(IgG1)

Sensorgrams




Fig 65. RBD-mut66-FC(IgG1)




Fig 66. RBD-mut67-FC(IgG1)




Fig 67. RBD-mut68-FC(IgG1)




Fig 68. RBD-mut69-FC(IgG1)

Sensorgrams




Fig 69. RBD-mut70-FC(IgG1)




Fig 70. RBD-mut71-FC(IgG1)




Fig 71. RBD-mut72-FC(IgG1)




Fig 72. RBD-mut73-FC(IgG1)

Sensorgrams




Fig 73. RBD-mut74-FC(IgG1)




Fig 74. RBD-mut75-FC(IgG1)




Fig 75. RBD-mut76-FC(IgG1)




Fig 76. RBD-mut77-FC(IgG1)

Sensorgrams




Fig 77. RBD-mut30-FC(IgG1)

Please note that the 'Fig Number' does not directly correspond to the 'mut number'. For example, Fig 77. corresponds to RBD-mut30-FC(IgG1).